Elektronische Lifterwaage

Electronic Lifter Scales

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Elektronische lift weegschaal

SR 282 0400

SR 232 1000

SA 630 0200

CE

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I: www.drivedevilbiss.fr

Table of contents

Foreword.....	19
Important safety recommendations.....	19
General data	20
Capacity and accuracy.....	20
Commissioning.....	21
Contents of the box	21
Attachment of the scale to the patient lift	23
Start-up of the scale.....	25
Function symbols and buttons.....	26
Operating anomalies.....	30
Cleaning.....	30
Technical data	31
Guarantee – Liabilities	31
Storage and transportation conditions	31
Handling and disposal of batteries.....	32
Manufacturer’s compliance statement.....	33

Foreword


This scale is a quality product that should be handled with care. Before using, please read the directions for use and follow the instructions provided. Retain this booklet.

Important safety recommendations

To ensure best performance and obtain accurate measurements, follow the usage and maintenance instructions set out within these directions for use.

- NEVER exert horizontal stress on the scale; there is a risk of breakage and irreversible malfunction of the weighing system.
- When moving, NEVER push the patient or the lifting-beam; use the patient lift handles.
- The scale and its fixings should be used EXCLUSIVELY on Drive DeVilbiss Healthcare France brand patient lifts.
- Comply with the maximum allowable weight for the patient lift. In addition, do NOT overload the scale. This could damage the sensors and would automatically invalidate the guarantee).
- Ensure that the patient is properly and safely installed within the sling.
- Ensure that the sling fixings are attached to the inside of the lifting-beam brackets and locked with safety clips.
- Do NOT leave any weight suspended on the scale for any length of time. This could damage the scale (reduction in the accuracy of measurements) and reduce the lifespan of the load cell.
- Before using the scale, ensure that the lower and upper fixings are correctly attached.
- When the display fades and is no longer very legible, think about replacing the batteries.
- Do NOT use the scale under environmental conditions (moisture, electronic interference ...) that could disrupt measurements.
- Do not attempt to open the scale. In case of problems with functioning, contact your retailer.
- Ensure that no cleaning fluid or water penetrates inside the scale.
- Dimensions (side in mm)

* Note: Dimensions are provided for information purposes only; they may change without prior notice.

General data

This scale is designed to be reliable and hard-wearing over time. It is easy to handle and simple to use.

Quality:

Design in line with a certified quality system – ISO 9001 standard.

Safety:

Stainless steel load cell with aluminium casing.

Product data:

20 mm LCD display, readable up to a distance of 10 m. 2 reliable attachment fixings.

Technology:

Technology including an electronic circuit and a reliable load cell, which ensures stability over time.

Power supply:

3 AA low-power batteries.

Handling:

Light and easy to transport.

Capacity and accuracy

Modell	Maximum Capacity (kg)	E	Division (n)
SR2321000 - HERMES	250 kg	0.1 kg	3,000
SR2820400 - SAMSOFT 175	175 kg		
SA6300200 - SAMSOFT and MINI SAMSOFT	150 kg		

Attention: Check the maximum allowable weight for the patient lift on which the scale is mounted.

Tare Range	100% (of the max. capacity of the scale)
Automatic re-set capability	$\pm 50\%$ (of the max. capacity of the scale)
Manual zeroing	$\pm 2\%$ (of the max. capacity of the scale)
Null seeking	0,5 E/s
Display read-out stability	≤ 10 seconds
Automatic stand-by	after 3 seconds (scale stable without button operation)
Automatic shut-down	after 3 minutes (scale stable without button operation)
Maximum overload	100% (of the maximum capacity of the scale) + 9 e
Maximum safe load	120% (of the maximum capacity of the scale)
Required usage temperature	- 10°C to + 40°C
Required usage humidity	= 90% at 20°C
Commissioning	

Insertion of batteries:

- Open the battery compartment located on the back by un-doing the screw for the cover.
- Insert three AA batteries, ensuring that the polarities (+/-) are correct or remove the insulator as necessary.
- Close the battery compartment, ensuring that it is closed properly.
- Replace the batteries if the indicator displays „LO“ = Batteries low.


Do not throw used batteries away – deposit them at an accredited collection point.

Contents of the box

- 3 AAA batteries
- 1 direction for use
- 1 scale
- Mounting set for Samssoft/Hermes

1 D10.2 int x D20 x L30 spacer (only provided for LP HERMES and which is mounted on the suspension axis at the boom end of the patient lift, so that the clevis rests on this brace)

1. Clevis x 1
2. Cotter pin x1
3. Threaded pin for a clevis x1
4. D10 int x D21 x L8 spacers x 2
5. Electronic scale (capacity 300 kg – accuracy 100 g) x 1
6. D10.3 int x D15 x L15 spacers x 2
7. Elastic ring x 4
8. Bracket x 2
9. Spacers (different depending on the patient lift version) x 2
10. Lifting-beam axle (different depending on the patient lift version) x 1
11. D10 L55 2-groove axle x 2

Attachment of the scale to the patient lift

Attachment of the scale should be carried out by an authorised and experienced technician; this operation is entirely the customer's responsibility and the greatest care needs to be taken with it. Drive DeVilbiss Healthcare France may not be held responsible in case of accident.

a) Preparation

Remove the lifting-beam from the patient lift boom, retaining all of the attachment elements (with the exception of the lifting-beam axle, which needs to be replaced by the specific pivot axle (10)).

b) Attachment of the clevis (1) to the top of the scale

Remove the safety pin (2), unscrew the clevis axle (3), place the clevis (1) either side of the upper fixing on the scale, then insert the clevis axle, placing the 2 washers (4) either side of the attachment and screw the clevis axle on fully. Place the safety pin (2) within the hole located at the end of the clevis axle, then bend the 2 ends of the pin.

c) Attachment of the lifting-beam pivot axle (10)

Use a flat-blade screwdriver to remove the two elastic rings from one side (7), then remove the bracket (8) and the two spacers (6 & 9). Insert the upper axle (11) and the spacer (6) into the lower fixing on the scale, replace the second spacer (6) from the other side of the attachment, insert into the upper hole on the bracket (8) and replace the elastic ring (7).

Insert the lifting-beam pivot axle (10) through the lifting-beam bearing, slide the lower axle (11) into the mounting hole on the lifting-beam axle, insert the second spacer (9), complete the action by sliding the axle into the lower hole on the second bracket (8), then re-attach the elastic ring.

ATTENTION:

Check that the 4 elastic rings are fully inserted into the grooves on the 2 axles (11); when this is the case, they should be able to turn freely around the axle.

For HERMES lifting-beam axles (10), the spacers (9) are chamfered; take care to position the chamfer against the lifting-beam pivot axle (10).

d) Attachment of the complete assembly to the end of the patient lift boom

Insert the lifting-beam axle into its housing at the end of the boom, then insert the clevis (1) and the original spacers (if there are any). Screw the nut back on, or replace the safety mechanism at the end of the boom axle.

Note: Depending on where the scale is mounted, it may be easier to start by attaching the scale to the end of the boom (step d), then finish by installing the lifting-beam under the scale (step c).

ATTENTION:

Before putting the patient lift back into service, ensure that the scale assembly and the lifting-beam are properly attached and secure.

The assemblies provided are intended exclusively for patient lifts manufactured by Drive DeVilbiss Healthcare France.

**Fitting of the sling**

Hang the sling on the 4 points of the lifting-beam. Ensure that the fixings are attached to the inside of the lifting-beam brackets and locked with safety clips

Start-up of the scale

Always ensure that the scale and the attachment fixings are correctly aligned vertically, with no constraints.

Press the  button to start up the scale. After a short test phase, the display should stabilise and indicate zero ,0.0' kg.

TARE function

The TARE button should be used to determine and deduct the weight of apparatus that will be supporting the patient (i.e. the weight of the lifting-beam, the sling ...).

- Firstly, weigh the equipment that will be carrying the patient (without the patient). When the weight is displayed, press the ,TARE' button to re-set the display to zero ,0.0'.

The weight of the equipment will be entered into the memory and will subsequently be deducted from the patient's weight during every weighing.

- Then, place the patient in the sling. Ensure that the 4 sling fixings are attached to the inside of the lifting-beam brackets and locked with safety clips.

- The patient's weight is displayed on the screen.

HOLD function

The ,HOLD' button enables the weighing result to be paused if there is any disturbance caused by the patient movements.

Function symbols and buttons

Start-up of the scale

Press  for 1 second to start the load test and calibration.

- ✓ The data ,STB' and ,ZERO' are displayed at the top left of the screen.

Shut-down of the scale

Press  for 1 second to turn the scale off.

- ✓ The screen displays the shut-down message ,off'.

TARE function

Press  to enter the tare value measured into the memory.

The symbol ,TARE' appears at the top of the screen and the value displayed returns to zero ,0.0'.

If the weight exceeds the maximum authorised load, or the value is negative, the ,TARE' function will not work.

If the weight is unstable, the ,TARE' function will become inactive.

As the tare weight is accounted for within the maximum load for the scale, this will be reduced accordingly.

Example: For a maximum load capacity of 300 kg, a tare of 100 kg will enable a weight of up to a maximum of 200 kg to be weighed.

Press  again to display the memorised tare value again.

The symbol ,TARE' disappears from the screen and the tare value is displayed.

HOLD function

Press  to pause in place the value displayed on the screen.

The symbol ,HOLD' appears and the value on the screen is paused in place.

Press  again to cancel the function.

The symbol ,HOLD' disappears from the screen and the value is no longer paused in place.

Units of measurement

Press  for 1 second to change the unit of measurement (3 possible choices: kg - kilogram, lb - pound or N - Newton).

If the ,TARE' function is active, the unit of measurement cannot be changed.

Re-set

Press  1 second to re-set the display to zero.
The symbol 'ZERO' appears at the top of the screen and the value displayed returns to zero '0.0'.

If a tare value is recorded, or if the scale is not stable or is held in place, a re-set cannot be carried out.

If a weight is suspended on the scale, a re-set cannot be carried out.

Adjustment of parameters

Press  and  simultaneously and hold these 2 buttons for 1 second to access the 'SETUP' Mode for adjustment.

Press  to access the 'off' automatic stand-by option and

press  to change the value.

The option may be set to 0 (deactivated), 5 (5 min), 15 (15 min), 30 (30 min) or 60 (60 min).

Press  to access the automatic Zero adjustment option and press  to change the value.

The option may be set to

- 0 (deactivated)
- 2 (2% of the max. capacity of the scale)
- 4 (4% of the max. capacity of the scale)
- 10 (10% of the max. capacity of the scale)
- 20 (20% of the max. capacity of the scale)
- 50 (50% of the max. capacity of the scale)

Press **TARE ZERO** to access the manual Zero adjustment option and press **HOLD UNIT** to change the value.

The option may be set to

- 0 (deactivated),
- 2 (2% of the max. capacity of the scale),
- 4 (4% of the max. capacity of the scale),
- 10 (10% of the max. capacity of the scale),
- 20 (20% of the max. capacity of the scale), or
- 50 (50% of the max. capacity of the scale).

Press **TARE ZERO** to access the screen back-lighting option and press **HOLD UNIT** to change the value.

The option may be set to ,off' (deactivated), 5 (5 sec), 15 (15 sec), 30 (30 sec), 60 (60 sec) or ,on' (continuous lighting).

Once the back-lighting has been set for a duration of 5, 15, 30 or 60 seconds, the count will start from the point at which the weight is stabilised.

When the back-lighting option is set to ,off', it is completely deactivated.

If the back-lighting option is set to ,on', the lighting will be continuously activated, without taking into account weight stabilisation.

Press **TARE ZERO** to exit the ,SETUP' Mode for adjustment.

Operating anomalies

Problem	Possible cause	Possible solution
Nothing is displayed on the screen after pressing the ON/OFF button	Batteries run down	Replace the batteries
	Batteries defective	Replace the batteries
	Incorrect insertion direction for batteries	Insert the batteries in the direction of the polarities
	ON/OFF button	Keep the button pressed down
The TARE or HOLD function does not work	TARE or HOLD button is defective	Clean the TARE or HOLD button
Unstable reading	Scale moving	Stabilise the load and the scale
The weight value measured seems to be incorrect	The scale was loaded before the test phase re-set to ,0.0' was completed The unit of measurement is not suitable	Do not load the scale. Wait for the display to indicate ,0.0' Check and select the desired unit of measurement

Cleaning

Clean with a damp cloth or with standard disinfection products.

Aggressive detergents should not be used.

Avoid allowing water or any other liquid to penetrate inside the scale.

Technical data

Accuracy	100 g
Dimensions	81 x 265 x 46 mm
Weight	Around 800 g (scale + fixings + batteries)
Power supply	3 x AA type alkaline batteries
Display	LCD-Display
Tare	Automatic Tare
Stand-by	Automatic stand-by

Guarantee – Liabilities

The product will be repaired or replaced if it is proved to have malfunctioned within a period of 1 year (one year) from the date of purchase (proof of purchase required), as a result of poor quality or material failure. Removable parts such as batteries are excluded from the guarantee. This guarantee does not cover abnormal use of the device, damage caused by falling, or that caused by incorrect use or inappropriate storage. In the case of equipment being returned, the scale should be packed in its original packaging; transportation costs will be payable by the purchaser. In case of any problem, you should contact the retailer from whom you purchased the equipment, whether this is because of a guarantee-related problem, the purchase of spare parts, or for repairs.

The guarantee will be void if the device is opened without authorisation.

Storage and transportation conditions

In the case of prolonged non-use, it is recommended that the batteries be removed and the scale be put away in its original packaging, or in any other packaging suitable for protecting it. Store everything in a clean and dry place.

It is advisable to retain the original packaging for the purposes of returning the equipment if need be, in order to avoid any damage during transportation. The guarantee will be void in the case of damage during transportation.

Handling of batteries

Do not deform or damage

Do not throw into the fire

Do not insert incorrectly

Do not open

Keep out of reach of children

Batteries do not belong in the household waste. As a consumer you are legally obliged to return used batteries. You can drop off your old batteries at public collection points in your community or wherever batteries of that type are sold.

NOTE:

These signs can be found on batteries:

Li = Battery contains lithium

Al = Battery contains alkali

Mn = Battery contains manganese

CR (Li); A (Al, Mn); AA (Al, Mn); AAA (Al, Mn)

Disposal of batteries

The symbol on the product or its packaging indicates that this product should not be treated as normal household waste but must be returned to a collection point for the recycling of electrical and electronic equipment. For more information, contact your local community, municipal waste management, or the shop where you purchased the product.

Manufacturer's compliance statement

This balance was manufactured according to harmonized European standards. It complies with the provisions of the EC Directives listed below:

- EMC Directive 2014/30 / EU on Electromagnetic Compatibility
- RoHS Directive 2011/65 / EC

in the applicable versions.

The explanation loses its validity, if at the Libra one not with us agreed amendment was made.

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